

Topspin Installation Guides

1. Get a copy of the software from us. Or go to the following website, register a free account, and then download the proper version for your system.
<https://www.bruker.com/service/support-upgrades/software-downloads/nmr/free-topspin-processing/free-topspin-download.html>
2. Install the software just like other software. Use the default settings through the installation. You will be asked to enter a superuser and password. Please use the administrator account for the PC.
3. Please check out the demonstration video (made by Dr. Keith Brown) for data processing using Topspin on our website. <https://sssc.usask.ca/nuclear-magnetic-resonance-laboratory.php>
4. Note the software has an expiration date. Please re-download and re-install the software after it is expired.